

Art. Nr.: MEG 125960

1.0 Zweck:

Zweck dieses Dokuments ist es, die Spezifikation der Sauerstoffbrille festzulegen.


2.0 Material

Farbe: Transparent / Maße Schlauch: Ø 5 mm x 2,1 m

3.0 Artikelbeschreibung

3.1 Verwendungszweck: Die Sauerstoffbrille ist für die Abgabe von Sauerstoffgas aus einer Quelle (Konzentrator, Tank, Regler usw.) an einen Patienten (über eine in die Nasenlöcher eingeführte Nasenkanüle) vorgesehen.

- 01 - Standardspitze
- 02 - Einlassrohr
- 03 - Schieber
- 04 - Y-Verbindungsstück
- 05 - Sauerstoff-Schlauch
- 06 - Sauerstoff-Anschluss



Materialübersicht



Artikel: Sauerstoffbrille, Standardspitze, Ø 5mm, 2,1 m Schlauch

Nr.	Bestandteil	Material	Nummer	Spezifikation
1	Standardspitze	PVC	1	Audlt
2	Einlassrohr	PVC	1	Ø 3 mm, 50 cm
3	Schieber	PP	1	/
4	Y-Verbindungsstück	PVC	1	/
5	Sauerstoff-Schlauch	PVC	1	Ø 5 mm, 2,1 m
6	Sauerstoff-Anschluss	PVC	1	/

Abmessungen der Nasenspitze

Maße	Höhe der Nasenspitze	Abstand zwischen Nasenspitze
Audlt	12-15 mm	8,4 mm<n<20,5 mm

Kanülen-Modell



1-125960-3-0001-2305

Specification of Nasal Cannula



Art. Nr.: MEG 125960

1.0 Purpose:

The purpose of this document is to define specification of the nasal cannula.


2.0 Material Specifications

Color: Transparent / Size Tube: Ø 5mm*2.1m

3.0 Product Descriptions:

3.1 Intended Use: Nasal Cannula is intended to deliver oxygen gas from a source (concentrator, tank, regulator, etc.) to a patient (via nasal nare inserted into nostrils).

- 01 - Standard Tip
- 02 - Inlet Tube
- 03 - Slide
- 04 - Y-Connector
- 05 - Oxygen Tube
- 06 - Oxygen Connector



Bill of Materials

Item: Nasal Cannula,Standard Tip, Ø 5mm, 2.1 m Tubing



No.	Part	Material	Number	Specification
1	Standard Tip	PVC	1	Audlt
2	Inlet Tube	PVC	1	Ø 3 mm, 50 cm
3	Slide	PP	1	/
4	Y-connector	PVC	1	/
5	Oxygen Tube	PVC	1	Ø 5 mm, 2.1 m
6	Oxygen Connector	PVC	1	/

Basic Dimensions of Nasal Tip

Size Height of nasal tip
Audlt 12-15 mm

Distance between nasal tip
8.4 mm<n<20.5 mm

Cannula Model



1-125960-3-0001-2305



megro GmbH & Co. KG
Am Schornacker 30,
46485 Wesel, Germany
www.megro.de



**Shanghai International
Holding Corp. GmbH (Europe)**
Eiffestrasse 80, 20537 Hamburg,
Germany



**Ningbo Shengyurui Medical
Appliances Co., Ltd.**
No. 138, BinHaisi Road, Hangzhou Bay New Zone,
Ningbo City, Zhejiang Province, 315336, P.R. China